



## **Intended Use**

The Para-Handy Bowling Wheelchair was designed to provide a lightweight means of mobility for users up to 110kg (17st). This wheelchair has been designed for the use on indoor and outdoor bowling greens.

The folding backrest and removable footrests make it easily transportable and for storage when not in use.

## **Technical Specification**

107cm	42"
66cm	26.5"
89cm	35"
52cm	20" - inc 2" cushion
40cm	16"
45cm	18"
40-50cm	16-20"
24" x 3"w	
3" Wide, t	o comply with (BWBA)
	66cm 89cm 52cm 40cm 45cm 40-50cm 24" x 3"w









## PARA-HANDY

## **BOWLING WHEELCHAIRS**

Approved by the BWBA, Bowls Scotland and the Scottish Indoor Bowling Association.

**Andy Semple** 

M: 07789028606

E: andrewgsemple@sky.com

Para-Handy Bowling Wheelchairs

www.para-handybowlingwheelchairs.com